

## APPENDIX IX PARAMETERS

CHEMICAL	TYPE	TEST CODE
Acenaphthlene	Semivolatile	CLR-A9-8270
Acenaphthylene	Semivolatile	CLR-A9-8270
Acetone	Volatile	CLR-A9-8260
Acetophenone	Semivolatile	CLR-A9-8270
Acetonitrile	Volatile	CLR-A9-8260
2-Acetylaminofluorene	Semivolatile	CLR-A9-8270
Acrolein	Volatile	CLR-A9-8260
Acrylonitrile	Volatile	CLR-A9-8260
Aldrin	Pesticide	CLR-8081
Allyl Chloride	Volatile	CLR-A9-8260
4-Aminobiphenyl	Semivolatile	CLR-A9-8270
Aniline	Semivolatile	CLR-A9-8270
Anthracene	Semivolatile	CLR-A9-8270
Antimony	Metal	A9-6010
Aramite	Semivolatile	CLR-A9-8270
Arsenic	Metal	A9-6010
Barium	Metal	A9-6010
Benzene	Volatile	CLR-A9-8260
Benzo[a]anthracene	Semivolatile	CLR-A9-8270
Benzo[b]fluoranthene	Semivolatile	CLR-A9-8270
Benzo[k]fluoranthene	Semivolatile	CLR-A9-8270
Benzo[ghi]perylene	Semivolatile	CLR-A9-8270
Benzo[a]pyrene	Semivolatile	CLR-A9-8270
Benzyl alcohol	Semivolatile	CLR-A9-8270
Beryllium (Total)	Metal	A9-6010
alpha-BHC	Pesticide	CLR-8081
beta-PHC	Pesticide	CLR-8081
delta-BHC	Pesticide	CLR-8081
gamma-BHC	Pesticide	CLR-8081
Bis(2-chloroethoxy)methane	Semivolatile	CLR-A9-8270
Bis(2-chloroethyl)ether	Semivolatile	CLR-A9-8270
Bis(2-chloro-1-methylethyl) ether AKA 2,2'-Oxybis(1-chloropropane)[Bis(2-chloroisopropyl)ether]	Semivolatile	CLR-A9-8270
Bis(2-ethylhexyl)phthalate	Semivolatile	CLR-A9-8270
Bromodichloromethane	Volatile	CLR-A9-8260
Bromoform	Volatile	CLR-A9-8260
4-Bromophenyl phenyl ether	Semivolatile	CLR-A9-8270
Butyl benzyl phthalate	Semivolatile	CLR-A9-8270
Cadmium (Total)	Metal	A9-6010
Carbon disulfide	Volatile	CLR-A9-8260
Carbon tetrachloride	Volatile	CLR-A9-8260
Chlordane	Pesticide	CLR-8081
p-Chloroaniline	Semivolatile	CLR-A9-8270
Chlorobenzene	Volatile	CLR-A9-8260

## APPENDIX IX PARAMETERS

CHEMICAL	TYPE	TESTCODE
Chlorobenzilate	Semivolatile	CLR-A9-8270
P-Chloro-m-cresol	Semivolatile	CLR-A9-8270
Chloroethane	Volatile	CLR-A9-8260
Chloroform	Volatile	CLR-A9-8260
2-Chloronaphthalene	Semivolatile	CLR-A9-8270
2-Chlorophenol	Semivolatile	CLR-A9-8270
4-Chlorophenyl phenyl ether	Semivolatile	CLR-A9-8270
Chloroprene	Semivolatile	CLR-A9-8270
Chromium (Total)	Metal	A9-6010
Chrysene	Semivolatile	CLR-A9-8270
Cobalt (Total)	Metal	A9-6010
Copper (Total)	Metal	A9-6010
m-Cresol	Semivolatile	CLR-A9-8270
o-Cresol	Semivolatile	CLR-A9-8270
p-Cresol	Semivolatile	CLR-A9-8270
Cyanide	Cyanide	CN-T-9014
2,4-D	Herbicide	A9-8151
4,4'-DDD	Pesticide	CLR-8081
4,4'-DDE	Pesticide	CLR-8081
4,4'-DDT	Pesticide	CLR-8081
Diallate	Semivolatile	CLR-A9-8270
Dibenz[a,h]anthracene	Semivolatile	CLR-A9-8270
Dibenzofuran	Semivolatile	CLR-A9-8270
Dibromochloromethane	Volatile	CLR-A9-8260
1,2-Dibromo-3-chloropropane	Volatile	CLR-A9-8260
1,2-Dibromomethane	Volatile	CLR-A9-8260
Di-n-butyl phthalate	Semivolatile	CLR-A9-8270
o-Dichlorobenzene	Semivolatile	CLR-A9-8270
m-Dichlorobenzene	Semivolatile	CLR-A9-8270
p-Dichlorobenzene	Semivolatile	CLR-A9-8270
3,3'-Dichlorobenzidine	Semivolatile	CLR-A9-8270
trans-1,4-Dichloro-2-butene	Volatile	CLR-A9-8260
Dichlorodifluoromethane	Volatile	CLR-A9-8260
1,1-Dichloroethane	Volatile	CLR-A9-8260
1,2-Dichloroethane	Volatile	CLR-A9-8260
1,1-Dichloroethylene	Volatile	CLR-A9-8260
trans-1,2-Dichloroethylene	Volatile	CLR-A9-8260
2,4-Dichlorophenol	Semivolatile	CLR-A9-8270
2,6-Dichlorophenol	Semivolatile	CLR-A9-8270
1,2-Dichloropropane	Volatile	CLR-A9-8260
cis-1,3-Dichloropropene	Volatile	CLR-A9-8260
trans-1,3-Dichloropropene	Volatile	CLR-A9-8260
Dieldrin	Pesticide	CLR-8081
Diethyl phthalate	Semivolatile	CLR-A9-8270
O,O-Diethyl phosphorothioate AKA Thionazin	Orthophosphated Pesticide	CLR-A9-8141
Dimethoate	Orthophosphated Pesticide	CLR-A9-8141
p-(Dimethylamino)azobenzene	Semivolatile	CLR-A9-8270
7,12-Dimethylbenz[a]anthracene	Semivolatile	CLR-A9-8270
3,3'-Dimethylbenzidine	Semivolatile	CLR-A9-8270
alpha,alpha-Dimethylphenethylamine	Semivolatile	CLR-A9-8270

## APPENDIX IX PARAMETERS

CHEMICAL	TYPE	TEST CODE
2,4-Dimethylphenol	Semivolatile	CLR-A9-8270
Dimethyl phthalate	Semivolatile	CLR-A9-8270
m-Dinitrobenzene	Semivolatile	CLR-A9-8270
4,6-Dinitro-o-cresol	Semivolatile	CLR-A9-8270
2,4-Dinitrophenol	Semivolatile	CLR-A9-8270
2,4-Dinitrotoluene	Semivolatile	CLR-A9-8270
2,6-Dinitrotoluene	Semivolatile	CLR-A9-8270
Dinoseb	Semivolatile	CLR-A9-8270
Di-n-octyl phthalate	Semivolatile	CLR-A9-8270
1,4-Dioxane	Semivolatile	CLR-A9-8270
Diphenylamine	Semivolatile	CLR-A9-8270
Disulfoton	Orthophosphated Pesticide	CLR-A9-8141
Endosulfan I	Pesticide	CLR-8081
Endosulfan II	Pesticide	CLR-8081
Endosulfan sulfate	Pesticide	CLR-8081
Endrin	Pesticide	CLR-8081
Endrin aldehyde	Pesticide	CLR-8081
Ethylbenzene	Volatile	CLR-A9-8260
Ethyl methacrylate	Volatile	CLR-A9-8260
Ethyl methanesulfonate	Semivolatile	CLR-A9-8270
Famphur	Orthophosphated Pesticide	CLR-A9-8141
Fluoranthene	Semivolatile	CLR-A9-8270
Fluorene	Semivolatile	CLR-A9-8270
Heptachlor	Pesticide	CLR-8081
Heptachlor epoxide	Pesticide	CLR-8081
Hexachlorobenzene	Semivolatile	CLR-A9-8270
Hexachlorobutadiene	Semivolatile	CLR-A9-8270
Hexachlorocyclopentadiene	Semivolatile	CLR-A9-8270
Hexachloroethane	Semivolatile	CLR-A9-8270
Hexachlorophene	Semivolatile	CLR-A9-8270
Hexachloropropene	Semivolatile	CLR-A9-8270
2-Hexanone	Volatile	CLR-A9-8260
Indeno(1,2,3cd)pyrene	Semivolatile	CLR-A9-8270
Isobutyl alcohol	Volatile	CLR-A9-8260
Isodrin	Semivolatile	CLR-A9-8270
Isophorone	Semivolatile	CLR-A9-8270
Isosafrole	Semivolatile	CLR-A9-8270
Kepone	Semivolatile	CLR-A9-8270
Lead (Total)	Metal	A9-6010
Mercury (Total)	Metal	HG-CV-7470
Methacrylonitrile	Volatile	CLR-A9-8260
Methapyrilene	Semivolatile	CLR-A9-8270
Methoxychlor	Pesticide	CLR-8081
Methyl bromide	Volatile	CLR-A9-8260
Methyl chloride	Volatile	CLR-A9-8260
3-Methylcholanthrene	Semivolatile	CLR-A9-8270
Methylene bromide	Volatile	CLR-A9-8260
Dibromomethane	Volatile	CLR-A9-8260
Methylene chloride	Volatile	CLR-A9-8260
Methyl ethyl ketone	Volatile	CLR-A9-8260
Methyl iodide	Volatile	CLR-A9-8260

## APPENDIX IX PARAMETERS

CHEMICAL	TYPE	TEST CODE
Methyl methacrylate	Volatile	CLR-A9-8260
Methyl methane sulfonate	Semivolatile	CLR-A9-8270
2-Methylnaphthalene	Semivolatile	CLR-A9-8270
Methyl parathion	Orthophosphated Pesticide	CLR-A9-8141
4-Methyl-2-pentanone	Volatile	CLR-A9-8260
Naphthalene	Semivolatile	CLR-A9-8270
1,4-Naphthoquinone	Semivolatile	CLR-A9-8270
1-Naphthylamine	Semivolatile	CLR-A9-8270
2-Naphthylamine	Semivolatile	CLR-A9-8270
Nickel (Total)	Metal	A9-6010
o-Nitroaniline	Semivolatile	CLR-A9-8270
m-Nitroaniline	Semivolatile	CLR-A9-8270
p-Nitroaniline	Semivolatile	CLR-A9-8270
Nitrobenzene	Semivolatile	CLR-A9-8270
o-Nitrophenol	Semivolatile	CLR-A9-8270
p-Nitrophenol	Semivolatile	CLR-A9-8270
4-Nitroquinoline 1-oxide	Semivolatile	CLR-A9-8270
N-Nitrosodi-n-butylamine	Semivolatile	CLR-A9-8270
N-Nitrosodiethylamine	Semivolatile	CLR-A9-8270
N-Nitrosodimethylamine	Semivolatile	CLR-A9-8270
N-Nitrosodiphenylamine	Semivolatile	CLR-A9-8270
N-Nitrosodipropylamine	Semivolatile	CLR-A9-8270
N-Nitrosomethylethylamine	Semivolatile	CLR-A9-8270
N-Nitrosomorpholine	Semivolatile	CLR-A9-8270
N-Nitrosopiperidine	Semivolatile	CLR-A9-8270
N-Nitrosopyrrolidine	Semivolatile	CLR-A9-8270
5-Nitro-o-toluidine	Semivolatile	CLR-A9-8270
Parathion AKA Ethyl Parathion	Orthophosphated Pesticide	CLR-A9-8141
PCBs	PCBs	CLR-8082
Polychlorinated dibenzo-p-dioxins	Dioxins	
Polychlorinated dibenzofurans	Furans	
Pentachlorobenzene	Semivolatile	CLR-A9-8270
Pentachloroethane	Volatile	CLR-A9-8260
Pentachloronitrobenzene	Semivolatile	CLR-A9-8270
Pentachlorophenol	Semivolatile	CLR-A9-8270
Phenacetin	Semivolatile	CLR-A9-8270
Phenanthrene	Semivolatile	CLR-A9-8270
Phenol	Semivolatile	CLR-A9-8270
p-Phenylenediamine	Semivolatile	CLR-A9-8270
Phorate	Orthophosphated Pesticide	CLR-A9-8141
2-Picoline	Semivolatile	CLR-A9-8270
Pronamide	Semivolatile	CLR-A9-8270
Propionitrile	Volatile	CLR-A9-8260
Pyrene	Semivolatile	CLR-A9-8270
Pyridine	Semivolatile	CLR-A9-8270
Safrole	Semivolatile	CLR-A9-8270
Selenium (Total)	Metal	A9-6010
Silver (Total)	Metal	A9-6010

## APPENDIX IX PARAMETERS

CHEMICAL	TYPE	TEST CODE
Silvex	Herbicide	A9-8151
Styrene	Volatile	CLR-A9-8260
Sulfide	Sulfide	H2S-9034
2,4,5-T	Herbicide	A9-8151
2,3,7,8-TCDD	Dioxin	
1,2,4,5-Tetrachlorobenzene	Semivolatile	CLR-A9-8270
1,1,1,2-Tetrachloroethane	Volatile	CLR-A9-8260
1,1,2,2-Tetrachloroethane	Volatile	CLR-A9-8260
Tetrachloroethylene	Volatile	CLR-A9-8260
2,3,4,6-Tetrachlorophenol	Semivolatile	CLR-A9-8270
Tetraethyl dithiopyrophosphate AKA Sulfotepp	Orthophosphated Pesticide	CLR-A9-8141
Thallium (Total)	Metal	A9-6010
Tin (Total)	Metal	A9-6010
Toluene	Volatile	CLR-A9-8260
o-Toluidine	Semivolatile	CLR-A9-8270
Toxaphene	Pesticide	CLR-8081
1,2,4-Trichlorobenzene	Semivolatile	CLR-A9-8270
1,1,1-Trichloroethane	Volatile	CLR-A9-8260
1,1,2-Trichloroethane	Volatile	CLR-A9-8260
Trichloroethylene	Volatile	CLR-A9-8260
Trichlorofluoromethane	Volatile	CLR-A9-8260
2,4,5-Trichlorophenol	Semivolatile	CLR-A9-8270
2,4,6-Trichlorophenol	Semivolatile	CLR-A9-8270
1,2,3-Trichloropropane	Volatile	CLR-A9-8260
O,O,O-Triethylphosphorothioate	Semivolatile	CLR-A9-8270
sym-Trinitrobenzene aka 1,3,5-Trinitrobenzene	Semivolatile	CLR-A9-8270
Vanadium (Total)	Metal	A9-6010
Vinyl acetate	Volatile	CLR-A9-8260
Vinyl chloride	Volatile	CLR-A9-8260
Xylene (Total)	Volatile	CLR-A9-8260
Zinc (Total)	Metal	A9-6010